

क्षेत्रीय आयुर्विज्ञान संस्थान, इंफाल: मणिपुर REGIONAL INSTITUTE OF MEDICAL SCIENCES, IMPHAL, MANIPUR 0385-2414750 mins edu.in (स्वास्थ्य और परिवार कल्याण मंत्रालय, भारत सरकार के अंतर्गत एक स्वायत्त संस्थान)

(An Autonomous Institute under the Ministry of Health & Family Welfare, Govt. of India)

Phone : 0385-2414720 website: www.rims.edu.in

NOTICE_INVITING OBJECTION

Subject:- Proprietary Certificate for Consumables & Reagents for Liquid Based Cytology, Brand: HOLOGIC, Model: ThinPrep 20000 Processor, for use in the Pathology Department, RIMS, Imphal - inviting comments/objections thereon.

This office is in the process for purchasing of Consumables & Reagents for Liquid Based Cytology machine for use in the Pathology Department, RIMS, Imphal, on Proprietary basis, from the firm M/s Hemogenomics Private Limited, Bangalore. The proposal submitted by M/s Hemogenomics Private Limited, Bangalore and PAC certificate are enclosed.

2. The above documents are being uploaded at RIMS website for an open invitation to all stakeholders to submit obligation/Objections or any comments, from the relevant manufacturers across India, regarding the proprietary and cost of nature of the equipment, within 15 days from the date of issue of this Notice, with reference to GEM Contract No.GEMC-511687722364518 dated 28.02.2024 for purchase of the Liquid Based Cytology machine for used in the Pathology Department, RIMS, Imphal. The comments should reach RIMS, Imphal, at email: rims.imphal@gov.in on or before 24th February, 2025 upto 4:00 P.M., failing which, it will be presumed that stakeholders have no objection to the Proprietary nature of the consumables and reagents items mentioned below:-

1. Table X - Thinprep Gym Test

Sl. no.	Particulars	Qty	Pack Size
i)	Preservcyt Int'l Solution Kit Asian (Gyn)	1 no.	250 vials/Box
ii)	Cervex Brush	1 no.	100 Pcs/Pack
iii)	Filter Thinprep Pap Test (Gyn)	1 no.	500 Filters/Box
iv)	Thinprep slide kit Gyn Asian	1 no.	500 Slides/Box

2. Table - Y Thinprep Non-Gyn Test

Sl. no.	Particulars	Qty	Pack Size
i)	Preservcyt Vial Non Gyn, Asia Pacific	1 no.	100 Vials Box
ii)	Cytolyt Quarts (4/Case) Asia Pacific	1 no.	4 Bottles/Pack
111)	Filter Thinprep Non-Gyn	1 no.	100 Filters/Box
iv)	Thinprep Microscope slides Non-Gyn	l no.	100 Slides/Box

3. Table Z - Thinprep strain reagents

Sl. no.	Particulars	Qty	Pack Size
i)	Thinprep Nuclear stain solution	1 no.	4 L/Pack
ii)	Thinprep Rinse solution	1 no.	4 L/Pack
iii)	Thinprep Orange-G-solution	1 no.	4 L/Pack
iv)	Thinprep EA stain solution	1 no.	4 L/Pack

1 no.

4 L/Pack

Encl: As stated above.

Sd/-(R.K. Mecolt Singh) Deputy Director (Admn.)

Copy to:
1. P.S. to the Director, RIMS, Imphal — for his kind information.

2. The HoD Pathology, RIMS, Imphal

3. The CAO/FA, RIMS, Imphal

4. The System Administrator, RIMS, Imphal — for uploading the tender notice with other relevant documents to RIMS website

Signed by Mecolt Rajkumar Singh Date: 14-02-2025 14:01:12

(R.K. Mecolt Singh) Deputy Director (Admn.)

Proprietary Article Certificate

(Procurement of the Goods under the provision of Rule 166 of GFR 2017 - Procurement base on Single Tender Enquiry)

The indenter hereby certify that

The indented goods/Software are manufactured/developed by

M/s HOLOGIC, Inc, localed at California 92121 USA

HEMOGENOMICS PRIVATE LIMITED, KARNATAKA

No other make or model/Software is acceptable for the following reasons: (ii)

Make: HOLOGIC

Model: THIN PREP 2000 and COMPASS STAINER

(iii) Supporting documents from supplier attached: Yes/No

Signature with date and designation of the indenting Officer

Name

D. SUSHMA KHURAIJAM

Signature

Designation of OFFICER

: Professor and HDD, Dept of Palhology : 23 Hold responsibility form

Dated

विकृति विद्यान विभाग

Department of Pathology क्षेत्रीय आगुविज्ञान संस्थान, इस्काल Regional Institute of Implied Medical Sciences, Imphal

LOLOGIC°

July 11, 2023

Proprietary Article Certificate

On behalf of Hologic, Inc., located at 250 Campus Drive, Marlborough, MA, USA 01.752, I hereby declare that the following are proprietary products, covered by intellectual property rights:

ThinPrep 2000 processor and ThinPrep 5000 processor utilizing a controlled membrane filtration technology that controls dispersion, collection, and transfer of diagnostic cells from the sample to the slide, in a closed system utilizing below reagents:

ThinPrep Gyn Kit:

PreservCyt Int'l Solution - Methanol-based reagent that serves as a transport, preservative, and antibacterial medium for gynaecologic samples.

Filter, ThinPrep Pap Test - Disposable plastic cylinder that is open at one end and has a filter membrane bonded onto the other end. The filter membrane has a flat, smooth, porous surface.

ThinPrep Slide kit Gyn Asian - It is a high-quality, pre-cleaned glass microscope slide.

ThinPrep Non-Gyn Kit:

PreservCyt 50 Pack Vial, Non Gyn, Asia Pacific - Methanol-based reagent that serves as a transport, preservative, and antibacterial medium for gynaecologic samples.

Cytolyt Quarts Asia Pacific (946 ml/Bottle)- is a water-based, buffered cell wash solution designed to lyse red blood cells, prevent protein precipitation, dissolve mucus, and preserve morphology of general cytology cellular samples.

Filter ThinPrep Non- Gyn Asian - Filters are single use consumable components intended to be used with a ThinPrep processor in the processing of fluid-based non-gynaecological cytologic specimens for examination. Blue in colour.

ThinPrep Microscope Slides Non-Gyn Asian - single use consumable components intended to be used with a ThinPrep processor in the processing of fluid-based non-gynaecological cytologic specimens for examination.

hinPrep Staining Kit:

ThinPrep Nuclear Stain Solution - is an aqueous solution that stains the nuclei of cells.

ThinPrep Rinse Solution - is an aqueous or alcohol-based solution that removes excess ThinPrep Nuclear Stain.

ThinPrep Orange G Solution - is an alcohol-based solution that serves to stain the cytoplasm of while maintaining nuclear staining and nuclear morphology.

ThinPrep Bluing Solution - is an aqueous, alkaline solution that causes cellular components stained with ThinPrep Nuclear Stain to darken to blue/violet.

ThinPrep EA Solution- is an alcohol-based solution that serves to stain the cytoplasm of cells on cytologic slide preparations.

Hologic, Inc. 250 Campus Drive, Marlborough, MA 01752 USA is the legal manufacturer and the owner of the intellectual property of the above-mentioned proprietary articles.

Signed on behalf of Hologic, Inc.

Name: Mark W. Irving

Title: Vice President and Secretary

Corporate Stamp:

- 2015 CERTIFIED

Mohan Cooperative Industrial Estate Matura Road, New Delhi – 110 044 Phone: +91-11-4132 5013, 4132 5014 Fax: +91-11-132 5015

17 June 2024

To,
The Director
Regional Institute of Medical Sciences
Lamphel Rd, Lamphelpat, Imphal
Manipur – 795004

Sub: Letter of Authority to Supply Consumables for Fully Automated High Throughput Liquid Based Cytology System

Dear Sir,

We, Hemogenomics Pvt. Ltd having office at 3rd floor, Site no. 26, Opp. Indian Petrol Bunk, New BEL Road, Near M.S. Ramaiah Hospital, Bangalore 560054 ("Hemogenomics") the authorized Distributors of Hologic Asia Pacific Limited who are established and reputable manufacturers of ThinPrep Products.

We, hereby authorize M/s Padma Medicare, having its office at RIMS Road, Nityananda Mandir, Imphal – 795001 (Manipur) with DL Number 13-X/DHS-2006 and GST Number 14BCNPS9876G1ZS to Supply Consumables for Fully Automated High Throughput Liquid Based Cytology System on our behalf.

We also hereby confirm that we would be responsible for the satisfactory execution of the supply.

With regards,

For Hemogenomics.PVT.LTD

Sangeeth G. Kini Vice President

Mobile: 09590302888

Hemogenomics Pvt. Ltd.

(9)